

URGENT MEDICAL DEVICE NOTIFICATION

BLUselect® Tracheostomy Tube Kits , BLUselect® Suctionaid® Tracheostomy Tube Kits, BLUgriggs® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with BLUselect® Tracheostomy Tube with or without Forceps, BLUperc® Dilation Procedural Tray with Single Stage Dilator Products, BLUperc® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with or without BLUselect® Tracheostomy Tube.

13 June 2024:


Dear Valued BLUSelect® Customers:

- Director of Biomedical Engineering
- Director of Nursing
- Director of Respiratory Care

Smiths Medical, Inc. is issuing this Urgent Field Safety Notice to notify you of a potential defect with the following BLUSelect®, BLUgriggs® and BLUperc® products listed in Attachment 1_Affected Product. This letter details the issue and the required steps for you to complete.

Issue:

Smiths Medical has identified the potential for a disconnection of the pilot balloon from the tracheostomy inflation line within specific lots of the BLUSelect®, BLUgriggs® and BLUperc® products because of a manufacturing defect. See the photo example of the issue below.



Potential Risk:

If the pilot balloon used to inflate the tracheostomy cuff becomes detached from the inflation line, the cuff pressure may not be maintained which can lead to inadequate ventilation of, and increased risk of aspiration to the patient. To date, Smiths Medical has received twelve (13) reports of serious injury and zero (0) deaths associated with this issue.

Affected Product:

See the Affected Product Codes Listed in the Attachment 1_Affected Product.

Smiths Medical Action:

Smiths Medical is sending this notification to all BLUSelect®, BLUgriggs® and BLUperc® customers who received product(s) from Smiths Medical listed in the Attachment 1_Affected Product. Smiths Medical will provide replacement product(s) and/or credit is available for affected customers. Please contact your local representative with your certificate of destruction to coordinate the replacement/credit.

Customer Required Actions:

When using the device, all instructions, including warning and cautions contained in the Instructions for Use documentation must be followed with heightened awareness. Please complete the following actions listed below:

1. Check all inventory locations within your institution for the specific lots of the affected product codes listed in the notification and discontinue use. Discard all affected products following your institution's process for discarding. If discarding is not immediately possible at your facility, then the product should be quarantined until disposal.
2. Share this notification with all potential users of the devices to ensure they are aware of this notification and proposed mitigations. If the devices are used at another location, please ensure this communication is delivered there.
3. Complete and return the attached Customer Response Form to smithsmedical8551@sedgwick.com ***within 10 days of receipt*** to acknowledge your understanding of this notification. Please contact your local representative for replacement product and/or a credit.
4. DISTRIBUTORS: If you have distributed potentially affected products to your customers, please immediately forward this notice to them. Request that they complete the response form and return it to smithsmedical8551@sedgwick.com

If you have any questions regarding product replacement and/or credit, please contact Smiths Medical Customer Service below:

Smiths Medical Contact	Contact Information	Areas of Support
Global Complaint Management	globalcomplaints@icumed.com 1-(866)-216-8806	To report adverse events or product complaints
Customer Service	Customerservice@icumed.com 1-(800)-0258-5361	Questions about product replacement and/or credit.

This notification is being performed with the knowledge of regulatory authorities, including the US Food and Drug Administration (FDA).

If you wish to contact the FDA regarding any adverse events or quality problems associated with this notice, use the following contact information.

- www.fda.gov/medwatch
- 1-(888)-INFO-FDA

Smiths Medical is committed to patient safety and is focused on providing exceptional product reliability and the highest level of customer satisfaction. Thank you for your prompt support on this important matter. We appreciate your cooperation.

Sincerely,



Andy Mathein
Vice President of Quality

Enclosures:

- Customer Response Form
- Attachment 1_Affected Product

101/817/080 BLUSELECT 8.0, CUFFED, FEN, W/WEDGE AND DECAN CAP

3933831 3962715 4062550 4068195 4230801 4294937

X

**101/596/ 080 BLUperc® Percutaneous Dilation Procedural Tray with 8.0mm
BLUselect® Suctionaid® Trach Tube**

3898770 3979602 4008832 3979615 4060531 4036527 4069871 4119536
4101970 4146207 4110837 4131437 4179935 4246523 4268459 4241328
4252236 4256036 4307794 4415596 4436549

X

101/595/090 BLUperc® Percutaneous Dilation Procedural Tray with 9.0mm BLUselect®

3898768 3944042 3950056 4008834 3987809 4036473 4036503 4008836
4060529 4098459 4168614 4110836 4177105 4183504 4291330 4294002
4324771 4330618 4263406 4265778 4344183 4423957

X

101/595/070 BLUperc® Percutaneous Dilation Procedural Tray with 7.0mm BLUselect®

3898766 4002180 4008828 4032081 4060523 4086681 4168564 4256056
4234424 4242669 4280697 4293992 4324235 4352941 4333064 4415609
4423960 4415588

X

101/892/080 BLUgriggs® Percutaneous Dilation Procedural Tray with forceps and

3943993 4165114 4226926

X

101/892/090 BLUgriggs® Percutaneous Dilation Procedural Tray with forceps and

39944020 4194294

X

101/892/070 BLUgriggs® Percutaneous Dilation Procedural Tray with forceps and

3950087

X

101/596/070 BLUperc® Percutaneous Dilation Procedural Tray with 7.0mm BLUselect®

3944058 3967805 4008830 4060560 4098464 4116821 4123503 4134674
4240514 4241327 4285358 4275910 4263399 4415589 4423966 4428303

X

101/596/090 BLUperc® Percutaneous Dilation Procedural Tray with 9.0mm BLUselect®

3950088 3944050 4008815 4032080 4060533 4119537 4098449 4101971
4123504 4261158 4280743 4285393 4270589 4263403 4301196 4319347
4423968 4326871 4436550 4119537 4098449 4101971 4123504 4261158
4280743 4285393 4270589 4263403 4301196 4319347 4423968 4326871
4436550 4447674

X

101/540/080 BLUgriggs Percutaneous Dilation Procedural Tray with forceps and

3943992 4258998

X

101/815/075 BLUSELECT 7.5, CUFFED, NON-FEN, W/WEDGE

X

101/815/080 BLUSELECT 8.0, CUFFED, NON-FEN, W/WEDGE, 1/EA

3905015 3940109 3942887 3965245 3965245 4021055
4058347 4029161 4100474 4118284 4187081 4169838
4190342 4202303 4204266 4206466 4192320 4204265
4274832 4209172 4294948 4336668 4336649 4336646

✓

101/815/070 BLUSELECT 7.0, CUFFED, NON-FEN, W/WEDGE, 1/EA

3905154 3949269 3965244 4039241 4032511 4091935 4115649 4121851
4187085 4131921 4187074 4196637 4219471 4214909 4220399 4223153
4277368 4260338 4306727 4306851 4292624 4294949 4306852 4332676
4332678 4332684 4336640 4369413

✓

101/595/080 BLUperc Percutaneous Dilation Procedural Tray with Single Stage Dilator

3898767 3950078 3950077 3987811 3991671 4008820 4008820 4032079
4032077 4060526 4072732 4086680 4098486 4110835 4131412 4148995
4128476 4163147 4179977 4265777 4256035 4272986 4283836 4296804
4324236 4326879 4330616 4339190 4259012 4341351 4425478

X

101/875/080 BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE, 1/EA

3905060 4038638 4071393 4091522 4106655 4083162 4088887 4124771
4220404 4187079 4187084 4202304 4279527 4294941 4375952

X

101/817/060 BLUSELECT 6.0, CUFFED, FEN, W/WEDGE AND DECAN CAP, 1/EA

3905152 4106973 4174794

X

101/815/100 BLUSELECT 10.0, CUFFED, NON-FEN, W/WEDGE

3933697 4029154 4148895 4112822 4161638 4021081 4187096 4169842
4192321 4284378 4206174

X

101/815/060 BLUSELECT 6.0, CUFFED, NON-FEN, W/WEDGE

3933722 3940999 3949223 4014730 4035539 4021067 4118283 4121883
4159150 4164109 4187097 4166710 4177295 4174793 4177293 4369411
4369412 4269270 4380967

✓

101/817/100 BLUSELECT 10.0, CUFFED, FEN, W/WEDGE AND DECAN CAP

3933799 4187080

X

101/815/060 BLUSELECT 6.0, CUFFED, NON-FEN, W/WEDGE

3905152 4106973 4174794

✓

101/815/090 BLUSELECT 9.0, CUFFED, NON-FEN, W/WEDGE

3933809 3949146 4011017 4026491 4076707 4021088 4159151 4169841
4185327 4178241 4200171 4169840 4204267 4200170 4279525 4212316
4216557 4260339 4294947 4334212 4336647 4336650

✓

3956297 4071378 4126686 4174782 4294945 4332682 4362485 4414625
4428012

101/815/075 BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE

3956339 4119879 4106981 4202297 4187077 4216554

X

101/815/085 BLUSELECT 8.5, CUFFED, NON-FEN, W/WEDGE

3956296 4090264 4071392 4209173 4332683 4369414

X

101/875/070 BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE

3965241 4040743 3997476 4040744 4040748 4083147 4090267 4106657
4106653 4130528 4219470 4244045 4244044

X

101/875/060 BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE

3973580 3994206 4077390 4077409 4106652 4097405 4121852 4174792
4185328 4228918 4294942 4306861

X

101/561/080 BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator

3887536 3906566 3906567 3906568 4288974 4285877 4291189 4348879
4348885 4308334 4348883 4348877 4386518

X

101/561/070 BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator

3928959 4005346 4217858 4234862 4266405 4348810

X

101/875/075 BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE

3965238, 4106984, 4230803 4318743, 4294936,

X

101/875/100 BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE

3905062

X

101/875/090 BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN, W/WEDGE

4034431 4068186 4067217 4106999 4106650 4115631 4118290 4187086
4206467 4216550 4220407

X

101/817/090 BLUSELECT 9.0, CUFFED, FEN, W/WEDGE AND DECAN CAP

3941008 4076692 4071380 4106651

X

101/817/070 BLUSELECT 7.0, CUFFED, FEN, W/WEDGE AND DECAN CAP

3986585 3956189 4107000 4118288 4187078 4130527 4208001

X

101/563/070 BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect®

3948836 4017238 4054722 4054801 4054723 4165746 4234758

X

101/561/090 BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator

3928993 3976736 4205592 4234863 4285761 4281078 4242949

X

101/540/090 BLUgriggs Percutaneous Dilation Procedural Tray with forceps and

4309778

X

101/563/080 BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect®

4348839

X